

Components

CATALOGUE

Accessories for the medical workplace

Precise solutions

With provita
components, you enable
the ideal design
of the medical
workplace.



Construction/Development

Developing medical products for the medical workplace in close cooperation with our partners and users.



Production at one location

Raw materials of the highest quality and processing in our own production facilities ensure efficient and quick production.



Products according to the Norm

We analyse, evaluate, produce and label our products in accordance with European directives and legal requirements.




Strong partnerships

We have been manufacturing quality products „Made in Germany“ for over 50 years. Trusting and long-term customer relationships are of great importance to us.

You are looking for an experienced partner that can help you realise your project?

Our project consultants with many years of experience are happy to advise and in cooperation with you, find specific solutions according to your wishes! Your planning will be also supported by drawings and 3D images.



Our sales team
Your contact person



Esteban Arco
e.arco@provita.de
+34 617 83 83 49



Christian Graf
c.graf@provita.de
+49 (0) 2193 5105-277



Michael Griep
m.griep@provita.de
+49 (0) 2193 5105-240



Micha Hilverkus
m.hilverkus@provita.de
+49 (0) 2193 5105-220



Lucica Krauss
l.krauss@provita.de
+49 (0) 2193 5105-210



Axel Meister
a.meister@provita.de
+49 (0) 2193 5105-230



Martin Vierkötter
m.vierkoetter@provita.de
+49 (0) 2193 5105-270



Henning vom Stein
h.vomstein@provita.de
+49 (0) 2193 5105-250



Marvin vom Stein
m.vomstein@provita.de
+49 (0) 2193 5105-260



Stephan Ziegler
s.ziegler@provita.de
+49 (0) 2193 5105-275



Your partner at provita: successful together as a team



»Through discussion and intensive exchange, improvements and new innovative products are created.«

Waldemar Worster,
with provita since 06/2004

» Quality as a lived matter of course. This applies to our products and also to the behaviour among each other «

Waldemar Mokrov,
with Provita since 08/2005



Wall rails

Medical rails for space-saving device management



IV-Poles

Professional infusion management on medical rails



Wall rail carts

Pivoting, mobile device carrier system for infusion and device management



Infusion and pump carrier systems

Practical and user-friendly infusion management



Lamps

Examination lamp
LED, halogen and twin version



Assembly & Accessories

Miscellaneous accessories for the extension of the wall rail system

Wall rails

The stable wall rail for flexible equipment

Flexible device management
on any wall, in any room.



High-quality
rail profile in
accordance to
an edge contour
standard





Stainless steel
(1.4301) for best
hygienic properties

MEDICAL RAIL
Load capacity
up to 100 kg
per metre

Material Stainless steel

Any length through connectors
(max. 4000 mm in one piece)

Load capacity 100.0 kg / m





MEDICAL RAIL
Easy installation with
screw-on discs

Material	Stainless steel
Cover and end caps	Plastic / black
Rail profile	25 × 10 × 1.5 mm
Load capacity	max. 100.0 kg (per metre)
Surface	sanded
Weight	approx. 1.5 kg (per metre)
Special feature	special and in-between lengths are also possible (with 1 grounding fixation at each rail)

Medical rail, please specify desired length Z2N27242



Plastic parts available
in white or blue colour
(without extra cost)

Optional cover
in stainless steel
(subject to sur-
charge) (Z2N06242)



Optional:
adaptable positioning
grounding box
(subject to surcharge)
(Z2N06132)



Connector (E2N4464-)



On customer request, we are also able to modify the distance between wall and rail.
Calculation of the number of bolts per rail = rail length / bolt space rounded to the nearest total number. Standard bolt space: 400 mm. Other bolt spaces on request (e.g: 500 mm or 600 mm – the load capacity is reduced with a bolt space of 500 or 600 mm). We reserve the right to provide rails larger than 3,000 mm with appropriate connectors.



For technical data
see provita.de

NEW

**MEDICAL RAIL WITH
CLIP PLASTIC COVER**
Easy to retrofit



huge
time saving of
approx. 1 min./
metre

Medical rail with clip plastic cover,
Please specify desired length Z2N27252



provita
COMPANY PRODUCTS BLOG NEWS SERVICE

Medical rail - planning tool

1. Enter general project information

Company name: Project name: Realisation time:

2. Add wall rails

Quantity: Length: Minimum length of a section: Minimum distance between the bolts: Minimum distance between bolts and rail end: Selection of covering caps: Rail number: (Rail number also entered separately by customer)

Add additional comments or special requests to this position

• Add position
• Delete last position
• Delete all positions

3. Submit configuration as working inquiry to provita

Customer number: Company name: First name: Last name: E-Mail address: Phone number: Street / House number: Postal code: City: Country:

☐ I have read the data protection information.


Send

easy to configure
online with our
planning tool:

[PROVITA.DE/RAILPLANNER](https://provita.de/railplanner)

IV-Poles

High-quality IV-Poles for every field of application



The matching wall rail IV-Poles
to your requirements.



Stepless height
adjustment



Stainless steel
(1.4301) for best
hygienic properties



Here you see two types of our
height adjustment (screw & single
handed system – the “original”)

INFUSION POLES Versatile & flexible

Material Stainless steel / Aluminium

Load capacity 2.0 kg / hook

Wall rail poles IV-Poles



+ different options

easy to
handle



Wall rail poles
IV-Poles



INFUSION POLE
With 2 wall
brackets

Material	Stainless steel
Tube diameter	25/18 mm
Length of outer tube	1,200 mm
Colour of covers	white
Angled	2 × 90°
Parallel offset	300 mm
Distance between wall and tube centre	120 mm
Height adjustment	Screw system
Bottle holder	Stainless steel 4 hooks
Regulating range	approx. 1,460 – 2,000 mm
Load capacity	2.0 kg / hook
Total weight	2.5 kg

I2002612

INFUSION POLE
Straight

Material	Stainless steel
Tube diameter	25/18 mm
Length of outer tube	1,000 mm
Height adjustment	Screw and Single-handed system
Bottle holder	Stainless steel 4 safety hooks
Regulating range	approx. 1,180 – 2,020 mm
Load capacity	2.0 kg / hook
Total weight	2.0 kg

Single-handed system

I1000422

Screw system

I2000422

INFUSION POLE
55° angled

Material	Stainless steel
Tube diameter	25/18 mm
Length of outer tube	1,000 mm
Angled	1 × 55°
Parallel offset	300 mm
Height adjustment	Screw system
Bottle bar	Stainless steel 4 safety hooks
Regulating range	approx. 1,250 – 1,800 mm
Load capacity	2.0 kg / hook
Total weight	2.9 kg

I2000572

INFUSION POLE
90° angled

Material	Stainless steel
Tube diameter	25/18 mm
Angled	2 × 90°
Parallel offset	300 mm
Height adjustment	Screw system
Bottle bar	Stainless steel 6 hooks
Load capacity	2.0 kg / hook
Total weight	3.0 kg

1,000 mm outer tube

Regulating range approx. 1,270 – 2,115 mm I2000552

1,500 mm outer tube

Regulating range approx. 1,780 – 2,400 mm I2000562

Wall rail poles
IV-Poles



INFUSION POLE
Stainless steel tubes
acc. to DIN EN 10217-7

Material	Stainless steel / Aluminium
Tube diameter	18/14 mm
Length of outer tube	800 mm
Height adjustment	Single-handed system
Bottle holder	Stainless steel 3 hooks – 2.0 kg / hook
Regulating range	approx. 1,150 mm – 1,700 mm
Total weight	1.9 kg
Equipped with	wall rail clamps, drip glass and holder, for mounting on medical rail 25 × 10 mm

Without wall support

I1000502

With wall support

I1000522

INFUSION POLE
With light-resistant
plastic parts

Material	Stainless steel / Aluminium
Tube diameter	25/18 mm
Length of outer tube	1,000 mm
Bottle holder	Stainless steel 3 hooks – 2.0 kg / hook
Regulating range	approx. 1,000 – 2,000 mm
Total weight	1.73 kg
Equipped with	wall rail clamp and wall support, for mounting on medical rail 25 × 10 mm

Single-handed system

I1002112

Screw system

I2000112

optional: for wall rail clamps and more see accessories



For technical data
see provita.de



INFUSION POLE
Angled

Material	Stainless steel / Aluminium
Tube diameter	25 / 18 mm
Length of outer tube	600 mm
Angled	2 × 45°
Parallel offset	approx. 350 mm
Bottle holder	Stainless steel 4 hooks – 2.0 kg / hook
Regulating range	approx. 1,180 – 1,510 mm
Total weight	1.45 kg
Equipped with	wall rail clamp and wall support, for mounting on medical rail 25 × 10 mm

Single-handed system

With additional screw height adjustment as anti-swivel protection	I1002212
---	----------

Screw system

I2002012

optional: for wall rail clamps and more see accessories

INFUSION POLE
With integrated pull-out stopper


Material	Stainless steel / Aluminium
Tube diameter	25 / 18 mm
Length of outer tube	800 mm
Bottle bar	Stainless steel 4 hooks – 2.0 kg / hook
Regulating range	approx. 1,000 – 1,600 mm
Total weight	2.03 kg
Equipped with	2 wall rail clamps, for mounting on medical rails 25 × 10 mm

Single-handed system

I1002712

Screw system

I2003112



INFUSION ROD
Safety hooks acc. to standards

Material	Stainless steel / Aluminium
Tube diameter	25 mm
Length	1,500 mm
Bottle bar	Stainless steel 4 hooks – 2.0 kg / hook
Total weight	1.5 kg
Equipped with	2 wall rail clamps, for mounting on medical rails 25 × 10 mm

I0001502

INFUSION POLE
With two extension clamps*

Material	Stainless steel / Aluminium
Tube diameter	38/25 mm
Length of outer tube	1,000 mm
Angled	2 × 90°
Parallel offset	350 mm
Height adjustment	Single-handed system with additional screw system as anti-swivel protection
Bottle bar	Stainless steel 4 safety hooks
Load capacity	2.0 kg / hook
Total weight	6.0 kg
Equipped with	2 wall rail clamps, for mounting on medical rails 25 × 10 mm*

I1002612

*Extensions clamps are available in other lenght as an option



ICU UNIT
Height adjustment
comfort

Material	Stainless steel / Aluminium
Height adjustment	Single-handed system
Bottle bar	Stainless steel 4 or 6 hooks – 2.0 kg / hook
Regulating range	370 mm
Total weight	approx. 2.1 kg

Equipped with wall rail clamp for mounting on medical rail or with mounting clamp for mounting on round tube 25 mm or 38 mm

With 4 hooks

For medical rail	I1009272
For round tube	I1009372

With 6 hooks

For medical rail	I1009292
For round tube	I1009392



INFUSION HOLDER
Space-saving
infusion supply

Material	Stainless steel / Aluminium
Bottle bar	Stainless steel 4 hooks – 2.0 kg / hook
Total weight	approx. 1.8 kg
Equipped with	wall rail clamp and wall support, for mounting on medical rail 25 × 10 mm

M0009004




implementing
individual
solutions



For technical data
see provita.de

Wall rail cart



Easy to change position with our mobile systems


For holding pumps, monitors
and other medical devices.

WALL RAIL CART Flexible positioning

Material Stainless steel / Aluminium


Load capacity 75.0 kg


Regulating range
1,100 – 1,900 mm



The wall rail systems can be
mounted in different ways
and different distances to
the wall.

Wall rail cart Wall rail cart





WALL RAIL CART A
With carrier tube

Carriage material	Aluminium
Carrier tube material	Stainless steel
Tube diameter	38 mm
Carriage colour	white
Length of outer tube	1,000 mm
System load capacity	75.0 kg
Equipped with	2 complete carriages and one carrier tube

Z2N0214A

Caution:
Wall rails must be mounted
exactly parallel to each other,
with clamping lever brake on
the upper carriage




Caution:
Please order stopper set
separately (Z2N0628A)

WALL RAIL CART B
With infusion pole

Carriage material	Aluminium
Carrier tube material	Stainless steel
Tube diameter	38 mm
Carriage colour	white
Length of outer tube	1,000 mm
System load capacity	75.0 kg
Total weight	approx. 4.2 kg
Equipped with	2 complete carriages and one infusion pole
Height adjustment	Single-handed system
Bottle bar	Stainless steel 4 hooks
Height adjustment	1,100 – 1,900 mm

Z2N0224A



WALL RAIL CART C
With angled infusion pole*

Carriage material	Aluminium
Carrier tube material	Stainless steel
Tube diameter	38 mm
Carriage colour	white
Length of outer tube	1,000 mm
System load capacity	75.0 kg
Total weight	approx. 4.2 kg
Equipped with	2 complete carriages and one angled infusion pole
Height adjustment	Single-handed system with additional screw system as anti-swivel protection
Bottle bar	Stainless steel 4 hooks
Height adjustment	1,325 – 1,850 mm

Z2N0225A


* to protect the devices from liquids



For technical data
see provita.de

Infusion and pump carrier systems

All devices within reach




For practical and
user-friendly infusion
& device management.



**MOUNTING
ON A WALL RAIL**
Horizontal




**MOUNTING
ON A ROUND TUBE**
Different diameters




**MOUNTING
ON THE WALL**

Material Stainless steel / Aluminium

Load capacity 2,0 kg / hook





**INFUSION
CARRIER SYSTEM**
Single

Material	Stainless steel / Aluminium
Regulating range	1,200 – 2,050 mm
Equipped with	2 pivoting support arms with universal clamp and a single-handed system infusion pole
Bottle holder	4 hooks
Support arm length	300 mm
Total weight	approx. 5.0 kg
I9900312	

**INFUSION
CARRIER SYSTEM**
Double

Material	Stainless steel / Aluminium
Regulating range	1,200 – 2,050 mm
Equipped with	2 infusion poles with single-handed system for mounting on a central tube Ø 38 mm
Bottle holder	4 hooks
Support arm length	300 mm
Total weight	approx. 7.9 kg
I9900412	


**EQUIPMENT
CARRIER SYSTEM**
Rotatable

Material	Stainless steel / Aluminium
Carrier tube length	1,000 mm
Equipped with	one support tube and two universal clamps
2 rotating arms	300 mm
Total weight	approx. 2.9 kg
Z2N6643S	

**DEVICE
CARRIER SYSTEM**
For vertical rails

Material	Stainless steel / Aluminium
Carrier tube length	700 mm
Equipped with	to mount on vertical 25 x 10 mm rails through height-adjustable brackets
Total weight	approx. 2.3 kg
Z2N6645A	

Can also be set
up as a vertical
support system



**EQUIPMENT
CARRIER TUBE**
Rotatable for wall
mounting

Material	Stainless steel / Aluminium
Carrier tube length	1,000 mm
Equipped with	one support tube and two wall brackets with two rotatable extension arms
Support arm length	200 mm
Total weight	approx. 3.0 kg
Z2N6647A	

DEVICE CARRIER TUBE
For wall mounting


Material	Stainless steel / Aluminium
Carrier tube length	1,000 mm
Equipped with	to mount on vertical 25 x 10 mm rails through height-adjustable brackets
Total weight	approx. 2.3 kg
Z2N6646A	



For technical data
see provita.de




Carrier systems
IWS and DUO-LOCK system



IWS UNIT
For infusion and syringe pumps
as well as infusion bottles

Material	Stainless steel / Aluminium
Regulating range	1,500 – 2,450 mm
Equipped with	one support tube Ø 33.7 mm, 2 horizontal swivel arms, 2 extension clamps for wall rail, 1 infusion pole with screw system, bottle bar with 6 hooks
Bottle bar	Stainless steel 4 hooks
Load capacity	2.0 kg per hook / 25.0 kg system
Total weight	approx. 9.0 kg
I9971702	



DUO-LOCK SYSTEM
Horizontal and
vertical adjustment

Material	Stainless steel / Aluminium
Regulating range	approx. 1,500 – 2,350 mm
Equipped with	2 Duo-Lock systems with quick hand fasteners and one single-handed system infusion pole
Bottle holder	Stainless steel 4 hooks
Load capacity	2.0 kg per hook / 10.0 kg overall system
I-DULO01	



For technical data
see provita.de

The all-round examination lamp



For technical data
see provita.de

For more lamps, please see our
lamp catalogue.



Made in
Germany

Lamps Series 1



**GOOSENECK
ARM LAMP**
Flexible positioning

Working radius approx. 800 mm

Total weight approx. 2.6 kg

Fixture pin Ø 18, L 35 mm

With spring balance joint

LED L100041A

Without spring balance joint

LED L100011A

Halogen L100005A

Twin version L100006A



**DOUBLE-JOINT
ARM LAMP**
Large projection

Working radius approx. 1,100 mm

Total weight approx. 2.9 kg

Fixture pin Ø 18, L 35 mm

With spring balance joint

LED L100081A

Halogen L100080A

Twin version L100082A

Without spring balance joint

Halogen L100061A

Twin version L100064A



**JOINT
ARM LAMP**
Simple and safe

Working radius approx. 800 mm

Total weight approx. 2.8 kg

Fixture pin Ø 18, L 35 mm

With spring balance joint

LED L100031A

Halogen L100030A

Twin version L100032A

Without spring balance joint

LED L100024A

Halogen L100020A

Twin version L100019A



RAIL CLAMP

To attach lamps
to 25 × 10 mm to
30 × 10 mm rails

Material Aluminium, ball-polished

Fixture entry Ø 18 mm


Weight approx. 0.2 kg

With threaded pin Z2N0643W



For technical data
see provita.de

Wall rails



Everything you need for mounting and the right accessories

For the completion and modification
of your wall rails.



EXTENSION CLAMP

For mounting on round tube
(Ø 25, Ø 33,7, Ø 38 mm)

Material Aluminium

Distance medical rail – round tube
100 mm

Colour white

Weight approx. 0.5 kg

For round tube Ø 38 mm Z2N09112

For round tube Ø 33.7 mm Z2N09122

For round tube Ø 25 mm Z2N09132



CLAMPS

For mounting lamps, poles and other
accessories, for mounting on medical
rail 25 × 10 mm

Material Aluminium, ball-polished

Weight approx. 0.4 kg

With grub screw, with fixture entry
Ø 18 mm Z2N0643W

With T-screw, with fixture entry
Ø 18 mm Z2N0647W

With T-screw, with fixture entry
Ø 25 mm Z2N0650W



WALL RAIL CLAMP KIT

With fixture entry for round tube
Ø 25 mm, for bridging unevenly deep
wall rails

Material Aluminium

Colour white

Weight approx. 0.8 kg

Special feature 1 adjustable
clamp (100 – 200 mm), 1 clamp with
fixed length

Z2N0618A



STOPPER KIT FOR WALL RAIL

For medical rail 25 × 10 mm

Material Aluminium

Colour white

For wall rail carts

Z2N0628A



CLAMP WITH QUICK-RELEASE LEVER

For mounting
on medical rail 25 × 10 mm

Material Aluminium, anodized

Weight approx. 0.1 kg

Z2N0932W



UNIVERSAL CLAMP

For medical rails and round tubes up
to Ø 40 mm, with four holes

Material Aluminium

Total weight approx. 0.3 kg

Z2N6638W



UNIVERSAL CLAMP WITH FIXTURE ENTRY BLOCK

For mounting on medical rails
25 × 10 mm or round tubes up to Ø 40 mm


Material Aluminium

Total weight approx. 0.45 kg

Fixture entry for round tube Ø 18 mm
Z2N6642S

Fixture entry for round tube Ø 25 mm
Z2N6641S

Holder for lamp pivot Z2N6639S



MEDICAL RAIL WITH UNIVERSAL CLAMP

For mounting infusion pumps and
other medical devices

Material Stainless steel / Aluminium

Profile 25 × 10 mm

Length 350 mm

Total weight approx. 0.85 kg

Z2N0684S



CARRIER TUBE WITH UNIVERSAL CLAMP

Universal holder for medical
equipment, for mounting on
horizontal medical rails

Material Stainless steel

Tube length 500 mm

Tube diameter 25 mm

Total weight approx. 0.9 kg

Z2N6644S



For technical data
see provita.de



MEDICAL RAIL
For round tube Ø 38 mm,
full length usable

Material	Stainless steel / Aluminium
Profile	25 × 10 mm
Length	558 mm
Weight	approx. 0,9 kg

Z2N06512



MEDICAL RAIL
680 mm
with 2 × clamp-sets for
Ø 38 mm round tube

Material	Stainless steel / Aluminium
Profile	25 × 10 mm
Centre distance	550 mm
Rail length	680 mm
Weight	1,3 kg

Z-DOCK05



MEDICAL RAIL
Retrofit mounting on round tubes
Ø 25 / 38 mm, usable length beside
the clamp 235 mm each

Material	Stainless steel
Profile	25 × 10 mm
Length	558 mm
Load capacity	approx. 8,0 kg
Weight	approx. 0,9 kg

For round tube Ø 25 mm Z2N70432
For round tube Ø 38 mm Z2N70332
For round tube Ø 33,7 mm Z2N70342



MEDICAL RAIL
double-sided 680 mm with
2 × clamp-sets for round tube Ø 38 mm

Material	Stainless steel / Aluminium
Profile	25 × 10 mm
Centre distance	550 mm
Rail length	680 mm
Weight	1,9 kg

Z-DOCK35



DOUBLE RAIL
For round tube Ø 25 mm

Material	Stainless steel / Aluminium
Weight	approx. 1,3 kg
optional	For round tube Ø 33,7 mm / Ø 38 mm
Rail length	558 mm

Z2N06932



DEVICE TRAY
For round tube Ø 38,
with edge protection

Material	Stainless steel / Aluminium
Size	300 × 250 mm
Load capacity	4,0 kg
Weight	approx. 2,75 kg

Z2N0754A



For technical data
see provita.de



SWIVEL ARM
With fixture pin
Ø 25 mm

Material	Aluminium
Centre distance	300 mm
Fixture entry	for Ø 25 mm round tube
Colour	white
Weight	approx. 1,0 kg

Z2N6650S



**SWIVEL ARM
FOR INFUSION POLE**
Clamp can be retrofitted

Material	Stainless steel / Aluminium
Centre distance	300 mm
Fixture entry	for Ø 25 mm round tube
Colour	white
Load capacity	25 kg

Z2N6653A



SUPPORT ARM
For holding
2 infusion poles Ø 25 mm

Material	Aluminium
Distance from tube entries to the centre point	300 mm
Tube support	Ø 25 mm, Ø 38 mm, Ø 25 mm
Colour	white
Weight	approx. 1,5 kg

Z2N6651S



SUPPORT ARM
For Ø 25 mm round tubes on both
sides

Material	Aluminium
Centre distance	250 mm
Fixture entry	for Ø 25 mm round tube
Colour	white
Weight	approx. 1,0 kg

Z2N6651A



**SWIVEL ARM
FOR INFUSION POLE**
Clamp can be retrofitted

Material	Stainless steel / Aluminium
Centre distance	300 mm
Fixture entry	for Ø 25 mm round tube
Colour	white
Load capacity	25 kg

Z2N6652A



PUMP CARRIER TUBE

Material	Stainless steel
Diameter	Ø 25 mm, Ø 38 mm

Ø 25 mm × 1000 mm	S0008172
Ø 25 mm × 1500 mm	S0008182
Ø 38 mm × 1000 mm	150411
Ø 38 mm × 1500 mm	291453

Assembly & Accessories
Wall rails



PUMP HOLDER

With wall rail clamp for mounting on medical rails 25 × 10 mm

Material	Stainless steel / Aluminium
Tube diameter	25 mm
Length	500 mm
Weight	ca. 0,95 kg

Z2N07162



WASTE BAG HOLDER

With wall rail clamp for mounting on medical rails 25 × 10 mm

Material	Stainless steel / Aluminium
Size	245 × 385 × 185 mm (W × H × D)
Weight	approx. 2.0 kg
Special feature	includes 3 waste bags

Z2N06202




BASKET WITH CLAMP

With wall rail clamp, for mounting on medical rails 25 × 10 mm, or with clamp-set to mount on Ø 38 mm round tube

Basket material	Steel, coated
Clamp material	Aluminium
Colour	light grey
Size	250 × 90 × 200 mm (W × D × H)
Weight	approx. 0.7 kg

for wall rail clamp
Z2S0665N

with clamp-set for round tube
Z-DOCK02



BASKET WITH 2 CLAMPS

With two wall rails clamps, for mounting on medical rails 25 × 10 mm

Material Korb	Steel, coated
Clamp material	Aluminium
Colour	light grey
Size	500 × 120 × 100 mm (W × D × H)
Weight	approx. 1.15 kg

Z2S0666N



TRAY

With two wall rail clamps, with edge protection

Material	Stainless steel / Aluminium
Tray size	300 × 250 mm
Load capacity	4.0 kg
Weight	approx. 2.75 kg

Z2N07332



GLOVE BOX HOLDER
3-FOLD FOR
MEDICAL RAIL

With wall rail clamp for mounting on medical rails 25 × 10 mm

Material	Stainless steel / Aluminium
Size outer dimensions	240 × 403 × 99 mm (W × D × H)
Size inner dimensions (per shelf)	240 × 132 × 96 mm (W × D × H)
Weight	approx. 2.5 kg

Z2N65522

Assembly & Accessories
Wall rails



GLOVE BOX HOLDER
2-FOLD
FOR MEDICAL RAIL

With wall rail clamp for mounting on medical rails 25 × 10 mm

Material	Aluminium, coated
Colour	silver
Size	280 × 270 × 100 mm (W × D × H)
Weight	approx. 1.3 kg

Z2N09353



GLOVE BOX HOLDER
3-FOLD
FOR MEDICAL RAIL

With wall rail clamp for mounting on medical rails 25 × 10 mm

Material	Aluminium, coated
Colour	silver
Size	390 × 260 × 100 mm (W × D × H)
Weight	approx. 1.5 kg

Z2N09363



CATHETER BASKET
FOR MEDICAL RAIL

Material	Stainless steel, electropolished
Mounting	Two hooks for medical rail
Weight	approx. 0.5 / 0.8 kg

150 × 100 × 280 mm (W × D × H)
Z2N08122

150 × 100 × 480 mm (W × D × H)
Z2N08132



CATHETER BOX
without separating screens

CATBOX 500 (height 500 mm)
Z2N08222

CATBOX 600 (height 600 mm)
Z2N08092

CATBOX 700 (height 700 mm)
Z2N08102

CATHETER BOX with two separating
screens (3-fold separation)

CATBOX 500 T (height 500 mm)
Z2N08272

CATBOX 600 T (height 600 mm)
Z2N08072

CATBOX 700 T (height 700 mm)
Z2N08082



CATBOX with window without screens

CATBOX 500 F (height 500 mm)
Z2N08232

CATBOX 600 F (height 600 mm)
Z2N08242

CATBOX 700 F (height 700 mm)
Z2N08252

CATHETER BOX with window and two
separating screens (3-fold separation)

CATBOX 500 FT (height 500 mm)
Z2N08262

CATBOX 600 FT (height 600 mm)
Z2N08052

CATBOX 700 FT (height 700 mm)
Z2N08062



For technical data
see provita.de



HOSE SUPPORT ARM

With gooseneck arm 360° rotatable

Material	Stainless steel/Plastic
Work radius	approx. 1190 mm
Fixture pin	Ø18 × L35mm
Weight	approx. 2,3 kg
Gooseneck arm	200 mm

Z2N06528




ARTICULATED ARM

With hose clamps (Ø approx. 23 mm),
without wall rail clamp

Material	Steel, chrome plated
Work radius	approx. 1,000 mm
Fixture pin	18 × L 35 mm
Weight	approx. 2.0 kg

Z2S06527




CABLE ARM

Double-extension to mount on rail 25 x 10 mm

Material	Stainless steel, electropolished
Centre distance	approx. 2 × 400 mm
Working radius	approx. 800 mm
Weight	1,5 kg

Z2N06832



STORAGE TRAY
FOR MEDICAL RAIL

Functional, space-saving foldable

Material	Stainless steel / Aluminium
Shelf dimensions	approx. 375 × 232 mm
Load capacity	approx. 10 kg
For medical rail	25 × 10 mm
Weight	approx. 2,6 kg

Z2N07602



SHELF

Without drawer

Material	Stainless steel
Shelf size	540 × 360 mm
Tragkraft	40.0 kg
Gewicht	ca. 6.0 kg

With two wall rail clamps and wall support Z2N09002

With two wall rail clamps and wall support,
with lateral rails Z2N09032

for wall rail cart

for round tube Ø 38 mm Z2N09042

for wall rail cart with lateral rails Z2N09162



SHELF

With 2 drawers

Tray material	Stainless steel
Drawer material	Steel, white coated
Shelf size	540 × 360 mm
Drawer size	420 × 310 mm
Drawer height	120 mm
Load capacity	40.0 kg



SHELF

With drawer

Tray material	Stainless steel
Drawer material	Steel, white coated
Shelf size	540 × 360 mm
Drawer size	420 × 310 mm
Drawer height	120 mm
Load capacity	40.0 kg

Weight approx. 10,1 kg, with 1 drawer Z2N09012

Weight approx. 10,8 kg, with 1 drawer
with lateral rails Z2N09182

for wall rail cart

Weight approx. 10,1 kg, for round tube Ø 38 mm,
with 1 drawer Z2N09152

Weight approx. 10,8 kg, for round tube Ø 38 mm
with 1 drawer and lateral rails Z2N09192

Weight approx. 14,7 kg, with 2 drawers Z2N09022

Weight approx. 15,4 kg, with 2 drawers
and lateral rails Z2N09212

for wall rail cart

Weight approx. 14,7 kg, with 2 drawers Z2N09172

Weight approx. 15,4 kg, for round tube Ø 38 mm
with 2 drawers and lateral rails Z2N09222





SHELF
With drawer


Tray material	Stainless steel
Drawer material	Steel, white coated
Shelf size	350 × 300 mm
Drawer height	120 mm
Load capacity	40.0 kg

With two wall rail clamps and wall support on request



OXYGEN BOTTLE HOLDER
for bottles of 2, 3 or 5 l with quick-release-lever

for bottles of 2 or 3 l	Z-DOCK28
for bottles of 5 l	Z-DOCK39



HOOK STRIP
To mount on Ø 38 mm round tube

Material	Aluminium / Stainless steel
Mounting	with doubled-clamp
Colour	white

2 hooks	Z2N3002A
2 hooks + 1 horizontal hook	Z2N3020A
3 hooks	Z2N3003A
3 hooks + 1 horizontal hook	Z2N3030A
4 hooks	Z2N3004A
5 hooks + 1 horizontal hook	Z2N3040A



OXYGEN BOTTLE HOLDER
height adjustable, for bottles of 2, 3 or 5 l
with quick-release-lever for wall rail 25 × 10
or with set of clamps for round tube Ø 38 mm

with quick-release-lever	Z-DOCK80
with set of clamps for round tube	Z-DOCK40



MEDICAL PLUGBOARD
“VARI MED”

Child proof (increased protection against
accidental contact with the plugs)

Housing material	Aluminium, anodized
Acc. to standards	DIN EN 60601-1-1

For mobile use
(cable size 1,5 mm²)

4-plug	Z2K06384
4-plug, dual fused, 3 m cable supply	Z2K07014
4-plug, dual fused, 5 m cable supply	Z2K07015
4-plug, spiral cable	Z2K07024
4-plug, dual fused, spiral cable	Z2K07034
6-plug, 3 m cable supply	Z2K06394
6-plug, dual fused, 5 m cable supply	Z2K07114
6-plug, dual fused, 3 m cable supply	Z2K07124



For vertical & horizontal mounting

variable
& robust



Create freedom of movement with infusion solutions from above.
You can find further information in our „Ceiling systems“ catalogue.

Instructions for Statics & Installation

Planning and building code requirements, as well as static tests, remain the responsibility of the buyer and must be verified before installation begins. Individual national and regional building codes must be considered. All products must be installed in accordance with commonly established technical rules.

Delivery does not include mounting parts for product installation as these depend on various static requirements. The operating instructions and safety instructions must be observed and complied with instructions for use must be kept in accordance with the legal regulations. Should the instructions for use not be enclosed by mistake, please contact us in any case.

Application of accessories

Always follow the instructions in the operating manuals of the product, when using accessories. If necessary, the maximum load of the accessories and the basic product must be reduced due to local conditions or the special configuration. In these cases there is no warranty granted by provita medical gmbh & co. kg. The compliance to the standards and regulations, especially the MPBetreibV (operator regulation) is the responsibility of the customer/user.

If you attach provita accessories to third party products or mount third party products to provita products, a new complete system may be created. The operator shall be then responsible for the evaluation and liability is then incumbent on the operator.

Installation of accessories

Generally when installing accessories and other parts, it is important to consider that:

- the mounting options are adapted to the product
- add-on components are firmly and securely attached
- the assembly of attachments does not create danger and potentially risks
- the balance point of the systems does not shift unfavourable (by offload parts)
- valid standards and guidelines are met (DIN ISO 60601, 3rd Edition, MPG, MPBetreibV)

Accessories for IV-Poles

We would like to emphasize that **not** all bottle holder/ - bar - models/ accessories are suitable for products with none weighted base constructions. For safety reasons, we would like to suggest to take advantage of the opportunity to get professional guidance by our competent sales team experts in regards to “tilting safety“ and the standard requirements of DIN EN 60601. Also changes done to the products jeopardise the approval in accordance with DIN EN 60601. Gladly we will provide consulting also in this matter.

Accessories for Wall and Ceiling Rails

When using accessories for Wall- and Ceiling Rail Systems, the specified maximum load capacity in the assembly instructions must be observed. The load capacity of the wall / ceiling must be calculated so that the input forces exerted by the rails can be absorbed. Relevant safety factors have to be considered. If applicable, static examinations are necessary.

General Product Information

provita products must only be used in accordance with their intended purpose, and in compliance with applicable statutory provisions, according to commonly accepted technical rules. They must be installed, operated and maintained in accordance with applicable work safety and accident prevention regulations.

They must be installed, operated, and maintained solely by persons with the **appropriate training, and adequate knowledge and experience.**

Before using any provita product, the user must ascertain its proper functionality and the conditions for effective application of the medical product. The user must also pay attention to the operating manual and other miscellaneous safety-related information and maintenance notes delivered with the product. A provita product operator must only charge persons, companies or institutions with the maintenance (maintenance, inspection, repair and reconditioning) of medical products, where the individuals concerned have the required technical knowledge, fulfill the necessary conditions, and have the necessary means to properly discharge such duties.

Upon maintenance or repair of medical products, the primary constructional and functional product features relevant to safety and operability must be inspected, in as far as they may have been affected by the maintenance activities performed.

Where safety-related inspections are not mandated, but at the same time are not expressly excluded, the operator shall, in accordance with generally accepted rules of the trade, perform or have such safety-related inspections performed at timely intervals, to enable the detection of any defects that may be expected in the light of past experience. The interval for such inspections shall, however, not exceed a period of two years.

An inspection log is to be created for these safety inspections, detailing the date of inspection and any evaluations made. The operator shall retain such logs at least until the

subsequent safety inspection. If a product is damaged or if there is any doubt as to its complete functionality, the product must not be used under any circumstances. In such specific cases, please feel free to contact us.

Disclaimer


provita medical's product liability shall be limited to the statutory liability requirements for the advertised characteristics of products contained in this catalogue. Any further liabilities – particularly for products from third parties listed in this catalogue – are excluded.

Trademarks

The products depicted in this catalogue may be under patent or subject to intellectual property rights. The reproduction of trademarked products constitutes a breach of the law, and may result in liability for damages. The information and illustrations in this catalogue do not release the user from the obligation to independently establish whether claims could arise from third party intellectual property rights.

Copyright provita medical gmbh & co. kg, 2020
Reproduction as well as duplication, also in extracts, are not permitted without permission of provita medical.

For more than 50 years, our family business has been one of the leading manufacturers in medical technology.





Our general terms and conditions apply. Errors, misprints and changes reserved.

© Copyright: provita medical gmbh & co. kg, november 2021

Reproduction as well as duplication, also in extracts, are not permitted without permission of provita medical.

**Browse also
in our other
product catalogues**



Version: 2021

**Auf der Huhfuhr 8
42929 Wermelskirchen
Germany**

**Phone +49 (0) 21 93 / 51 05 - 0
Fax +49 (0) 21 93 / 51 05 - 269**

**info@provita.de
www.provita.de**

